



ABSTRACT


Altratech introduce a novel method of molecular detection. This method utilizes synthetic peptide nucleic acid (PNA) probes for NA capture, antigenic peptides for antibody capture, superparamagnetic beads, electromagnetics, microfluidics and CMOS silicon chip technology. Optical detection is replaced with electronic detection. This novel assay enables simultaneous genetic and serology testing.

Results are presented for detection of HIV-1 & HIV-2 antibodies from HIV-positive patient plasma; and detection of SARS-CoV-2 RNA virus from human saliva. Altratech's technology, when encapsulated in inexpensive single use cartridges, will enable complex molecular testing to be undertaken at the point of need^{1,2}.

PROPRIETARY MATERIALS


Digital CMOS Sensor Chip³:

Single paramagnetic beads are specifically tethered by the target to our sensor & quantitatively detected by employing fringe-field sensing:



PNA Probes (for NA capture):

Altratech have designed PNA probes for specific target capture directly from samples and have co-developed chiral PNAs synthesized in Fmoc chemistry with the US National Institute of Health⁴. PNA probes are known for their excellent specificity, e.g. Single-Base-Mismatch specificity:




Antigenic Peptides (for antibody capture):


Synthesized with the same Fmoc chemistry as PNA.

Superparamagnetic Beads:

PNA-coated superparamagnetic beads with captured Target are magnetically removed from the sample⁵.




PROTOTYPE & HARDWARE



2025 FULLY INTEGRATED DEVICE

- Full integration of sample-input with detection is underway in 2025, in conjunction with Cambridge Design Partners (UK)
- Clinical validation with 140 HIV patient samples is planned at St. Cecilio Hospital, Granada, Spain.




2026 HIGH VOLUME PRODUCT

The final product will be a small form-factor, low-cost cartridge and wireless reader unit (Fig 6). This miniaturisation is made possible through the solid-state nature of the device and by replacing Optical Detection with our CMOS enabled Electronic Detection:



ASSAY METHOD⁵

There are two assay variants: (1) Indirect and (2) Direct.




The indirect approach facilitates rapid progress towards improved assay sensitivity, through electromagnetic mixing and hardware improvements. As upstream advances in Target Enrichment are implemented, the direct method (no replication step) will supplant the indirect approach when the sensitivity levels are comparable.

INDIRECT: NUCLEIC ACID DATA

- **Current LOD** 1E4 copies/ml using SARS-COV-2.
- **Targeting** 1E3 copies/ml by end of Q1 2025. This is the WHO clinically relevant 1000cp/mL level for portable non-laboratory RNA viral-load testing¹.
- **Capable** of achieving 1E2 copies/ml sensitivity with further engineering refinement.

DIRECT: NUCLEIC ACID DATA


- **Current LOD** 1E6 copies/mL using synthetic SARS-COV-2.
- **Targeting** 1E4 cp/mL by end of Q3 2025.
- **Capable** of achieving 100 cp/mL sensitivity with further engineering and target enrichment development.



DIRECT: ANTIBODY DATA

The key principle of our technology, **Detection by Proxy⁵**, enables this hardware and **direct** assay to also **detect and quantify protein targets**:

Fig 8: Capacitance-vs-HIV1 & HIV2 immuno-assay **Antibodies** (n=2).
→HIV 1 & 2 antibodies detected in HIV-positive-patient plasma (NIBSC).



Here we demonstrate the replacement of the enzyme in a commercially-available HIV antibody test with our superparamagnetic PNA coated reporter beads. This allows specific detection by complementary PNA's on the sensor chip.

DISCUSSION

Analytical specificity:

PNA's neutral charge enables them to bind targets in raw samples, with **excellent specificity** (Fig.2). This allows for sequence-specific extraction and eliminates centrifuging and complicated sample-prep steps from this assay. Many different infectious diseases can be detected simply by changing the PNA probe sequences.

Analytical sensitivity & LOD:

The *indirect* method in Q1'2025 is facilitating rapid progress to improved sensitivity. The *direct* method will enable further improvements and cost-reduction, to target compliance with the challenging demands of the WHO REASSURED criteria⁶, most especially in providing a commercially viable diagnostic tool for the developing world.

CONCLUSION

- We have introduced a first-of-kind multi-analyte combo DNA/RNA/Antibody/Antigen diagnostic.
- This technology can be used for both Serology and Nucleic Acid detection and quantification. Furthermore, it offers capability to multiplex multiple targets in the same assay.
- Its stand-alone portability makes it ideal for point of care triage and point of need testing.
- Applications for this technology are in molecular diagnostics, companion diagnostics and vaccine development.

REFERENCES

1. HIV diagnostic rec, WHO, ISBN 978-92-4-151621-1
2. www.msfacecess.org/time-for-5
3. "Next Generation Molecular Detection with a Capacitive Sensor", T. Cummins, B. O'Farrell, (doi:10.1007/978-3-031-28912-5_6)
4. Cyclopentane Peptide Nucleic Acid (PNA)", D. Appella, B.O'Farrell, K.Oshaben et al, Biopolymers, (doi.org/10.1002/bip.23481, Biopolymers 2021)
5. Altratech patents US11274291/10738348/11796498, US11459601 and EU/CN/JP equivalent granted patents. (Altratech Ltd owns an extensive patent portfolio protecting all aspects of this product & technology).
6. https://www.nature.com/articles/s41564-018-0295-3

ACKNOWLEDGEMENTS

- Tech-transfer of PNA synthesis from NIH to Altratech & SARS-CoV-2 capture from saliva was funded by **BARDA contract 75A50122C00067**.
- End-to-end assay proof-of-principle was funded by **EU-H2020 'Easiviral' grant 768090** & Enterprise Ireland.
- The 2025 'Triple-H' HIV-HCV-HBV production is being funded by the European-Innovation-Council **EIC-Accelerator 'UbiHealth'** grant (101188322) and including EU public-private matched equity investment.